



OPTIMUM GROWTH

TIPS & TRICKS FOR CHO & HEK



Thomson Cell Culture Solution

Thomson Instrument Company was founded in 1970 to service the needs of the growing life sciences industry. It remains a family run business with a long history of product innovation. In addition, it should be noted that all Thomson products are designed and manufactured in the USA. Thomson has worked with its customers to design products that outperform the competition by introducing unique design features that add functionality to an otherwise static product market. In particular, Thomson has reinvented shake flasks that are used for culturing *E.coli* and mammalian cells for recombinant DNA and protein production. These features are

reviewed in Dr. Bruecher's article, "Evolution of Shake Flask Technology. Novel Product Introductions Offer Advantages by Increasing DNA and Protein Production."¹

The advantages of Thomson Optimum Growth™ Flasks (patented) can be attributed to a fairly simple principal of good mixing and high oxygen transfer rate helping cells grow to higher densities. Increased aeration is achieved through a patented baffle design along with a large diameter neck and gas permeable cap.

The innovation in Thomson's Optimum Growth™ product line does not end with a

single flask. Thomson has worked with its customers to reduce the footprint of its flasks in such a way to maximize space in the most commonly used cell culture incubators.

In the tradition of evolving to fit customer needs, Thomson has introduced a line of companion products that work along with the Optimum Growth™ Flasks. These include caps, lids & ports that allow the sterile transfer of fluid into and out of the flasks. Current cell culture processes often involve longer culture durations that require both sampling over time and batch feeds. Thomson products support these existing workflows and help to streamline the process wherever possible.



Optimum Growth™ Flasks








Optimum Growth™
Transfer Caps



Optimum Growth™
Special Flasks

¹Bruecher. September 15, 2015. Novel Product Introductions Offer Advantages by Increasing DNA and Protein Production. Genetic Engineering & Biotechnology News. 35 (16). <http://www.genengnews.com/issue/toc/255>

Page Number	Product	Cell Growth	Aseptic Gravity Transfer	Aseptic Pump Transfer	Sampling	Feeding	0.2µm PTFE Vent Cap	Sterile
4	 Optimum Growth™ Flasks Available in 125mL - 5L for insect & mammalian cell growth.	●					●	●
14	 Sampling Flasks Available in 125mL, 250mL, 500mL & 5L Optimum Growth™ Flask sizes. Provides simple aseptic way to sample cells in the shaker.	●			●		●	●
17	 Multiported Optimum Growth™ Flasks Available in 1.6L & 5L with ports for aseptic transfer, feeding & sampling during cultivation.	●		●	●	●	●	●
19	 Inversion Transfer Caps Gravity transfer from 1.6L, 2.8L & 5L Optimum Growth™ Flasks to flasks, bags or bioreactors. Removing the peristaltic pump from the process reduces stress on cells during transfer.		●				●	●
22	 Bidirectional Transfer Caps Sterile bidirectional transfer utilizing a peristaltic pump from 1.6L, 2.8L & 5L Optimum Growth™ Flasks to flasks, bags or bioreactors or to feed/seed flasks during cultivation.			●		●	●	●
24	 Rapid Clear® 3000 Cap An innovative new technique that allows high speed clarification of cellular material	●		●				●



Key Features

- Baffles designed for High Aeration & Low Shear
- Same Footprint as Comparable Fernbach Flask
- Less Foaming than Disposable Fernbach
- Transfer Cap connects directly to bags & bioreactors with quick connect or tube fusing
- .2µm Vented Cap
- Individually Packaged and Sterilized

OptimumGrowth™


The advantages of Thomson patented Optimum Growth™ Flasks can be attributed to a fairly simple principal of good mixing and high oxygen transfer rate helping cells grow to higher densities. Increased aeration is achieved through a patented baffle design along with a large diameter neck and gas permeable cap. The innovation in Thomson's Optimum Growth™ product line does not end with a single flask. Thomson has worked with its customers to reduce the footprint of its flasks in such a way to maximize space in the most commonly used cell culture incubators. In the tradition of evolving to fit customer needs, Thomson has introduced a line of companion products that work along with the Optimum Growth™ Flasks. These include caps, lids & ports that allow the sterile transfer of fluid into and out of the flasks. Current cell culture processes often involve longer culture durations that require both sampling over time and batch feeds. Thomson products support these existing workflows and help to streamline the process wherever possible.



**Optimum Growth™ Flasks
have 0.2µm PTFE Vent Caps**

Optimum Growth™ Flask Specifications

Flask Size	125mL	250mL	500mL	1.6L	2.8L	5L
Part #	931110	931111	931112	931113	931114	931116
Image						
Description	Optimum Growth™ 125mL Flask w/ Vent Cap Sterile	Optimum Growth™ 250mL Flask w/ Vent Cap Sterile	Optimum Growth™ 500mL Flask w/ Vent Cap Sterile	Optimum Growth™ 1.6L Flask w/ Vent Cap Sterile	Optimum Growth™ 2.8L Flask w/ Vent Cap Sterile	Optimum Growth™ 5L Flask w/ Vent Cap Sterile
Qty/Case	50	50	25	12	6	4




Thomson Optimum Growth™ Flasks are designed in a way that protein production scales consistently across all sizes, unlike any other shake flasks on the market. The Optimum Growth™ Flask Family has a unique patented design that allows one shake speed for all sizes. Multiple flask sizes on one shaker platform eliminates the need for additional shakers at different speeds.

Consider the Health of Your Cells:

Complex Cellular Systems

Are your cultures healthy? Cellular stress can be triggered by environmental changes including oxygen/CO2, pH, and acetate levels. These stresses will lower the viability and protein yield of cultures. The Optimum Growth™ Flasks patented design improve aeration improving both quality and quantity of the cell culture.



Comparison of 293 & CHO Cell Line Traits for Transient Expression

Transient expression of HEK-293 and CHO-K1 cell lines are frequently used for antibody production. Both cell lines are easily adaptable to suspension culture, and grow to high cell densities in the Optimum Growth™ Flasks.

Characteristic	293	CHO
Tolerance to culture stress	Best	
Ability to shift from growth to expression		Best
Speed of expression	Best	
Duration of expression		Best
Stability of transfection complexes	Best	
Lower amount of DNA needed		Best
Flexibility of protocol		Best
Benefit of batch feeds		Best

Excellent Growth with Space Saving Capability

Compared to competing flask systems, Thomson's Optimum Growth™ line of flasks have higher working volumes while keeping the same footprint. Using smaller flask sizes with similar working volumes to larger competing flasks allows a higher total volume per shaker.



18 x 1.6L Optimum Growth™ Flasks
Total Volume 16.2L/shaker



12 x Optimum Growth™ 2.8L Flasks
Total Volume 16.8L/Shaker



8 Position Carrier available
for 125mL & 250mL Flasks


Fill Volumes & Shake Speeds

For all tables | 1" = 25mm | 2" = 50mm

CHO Stable Cells, CHO Transient, HEK 293 Transient , Hybridoma					
Best Fill Volume	Working Volume	Flask Size	Best Speed *RPM in 1"/2"	Speed Range 1"	Speed Range 2"
63mL	24-75mL	125mL	150 / 110	120-150	90-110
150mL	100-150mL	250mL	150 / 110	120-150	90-110
250mL	175-250mL	500mL	150 / 110	120-150	90-110
900mL	0.4-1.1L	1.6L	150 / 110	120-150	90-110
1.4L	0.9-1.6L	2.8L	150 / 110	120-150	90-110
2.5L	1.7-3.2L	5L	120 / 90	120-140	90-110

Insect Cells					
Best Fill Volume	Working Volume	Flask Size	Best Speed *RPM in 1"/2"	Speed Range 1"	Speed Range 2"
63mL	24-75mL	125mL	150 / 110	120-150	90-110
150 mL	100-150mL	250mL	150 / 110	120-150	90-110
250mL	175-250mL	500mL	150 / 110	120-150	90-110
900mL	0.4-1.1L	1.6L	150 / 110	120-150	90-110
1.4L	0.9-1.6L	2.8L	150 / 110	120-150	90-110
2.5L	1.7-3.2L	5L	135 / 95	120-140	90-110

Microbes/ <i>E.coli</i>			
Best Fill Volume	Working Volume	Flask Size	Best Speed *RPM in 1"/2"
63mL	24-75mL	125mL	250 / 150
125mL	100-150mL	250mL	250 / 150
250mL	175-250mL	500mL	250 / 150
900mL	0.4-1.1L	1.6L	250 / 150
1.4L	0.9-1.6L	2.8L	250 / 150
2.5L	1.7-3.2L	5L	250 / 150



Optimum Growth™ Flask FAQs

What have people done successfully to change vessels from Spinner flasks & Roller bottles to Optimum Growth™ Flasks (patented)?

Cells adapted to Spinner Flasks and Roller Bottles can be easily transitioned to Optimum Growth™ Flasks. Adjusting existing cultures from different formats to Optimum Growth™ Flasks requires reducing the volume and shake speeds of the first 1-2 passages*. The addition of up to 1% of surfactant** to the media may be needed due to spinner flasks and roller bottles having lower shear than shake flasks. Once the cells have adjusted to the shake flasks, recommended speeds will work well.

* See chart with fill volumes and shake speeds.

** ThermoFisher Pluronic, p/n 24040032 or MilliporeSigma Simethicone, p/n 59920C

Why do Optimum Growth™ Flasks work better than other disposable flasks (non-baffled or baffled) for mammalian cell lines (CHO, HEK293, etc.) & insect cell lines (SF-9, SF-21, High Fives, Trichoplusia ni)?

Optimum Growth™ Flasks are patented shake flasks designed for high aeration and low shear. Optimum Growth™ Flasks achieve high aeration due to a unique baffle design that has been optimized for mammalian and insect cell lines. They provide enhanced gas exchange with low shear mixing, which can increase yields significantly when combined with both nutrient enriched media and proper pH balance.

What clamps and shakers work best with the Optimum Growth™ Flasks?

Optimum Growth™ Flasks are designed to shake in 1" or 2" orbit shakers. Sticky tape or rug gripper pad is recommended for under 170rpm. Our 125mL, 250mL and 500mL flasks will work with standard shake flask clamps.

Thomson Flask Clamp Compatibility Chart

Flask Size	Eppendorf®	INFORS HT®	Kuhner®	Fisher Scientific®	VWR®
5L	ACE-5000S	12207 ATM 3000	SM313000F	50-901-18155	57019-696
2.8L	ACE-2000S	ATK 2000	SM312800T	14-281-27	57019-686
1.6L	N/A	31272	SM311600	N/A	N/A
500mL	M1190-900	12204 ATM 500	SM310500	50-901-18152	N/A
250mL	M1190-9002	12203 ATM 250	SM310250	50-901-18151	N/A
125mL	M1190-9001	12202 ATM 100	SM310125	50-901-18150	N/A

Are the Optimum Growth™ Flasks single use?

Yes, the Optimum Growth™ Flasks are designed for single use and are not autoclavable. They are competitively priced compared to disposable bioreactors or shake flasks from other manufacturers.

What are the Transfer Caps that go along with the Optimum Growth™ Flasks?

Inversion & Bidirectional Optimum Growth™ Transfer Caps (patented) allow for a quick stress free cell transfer between flask and downstream vessel (Optimum Growth™ Flasks, cell culture bags, bioreactors, etc.). Inversion Transfer Caps use the power of gravity to facilitate transfer, thus maintaining higher culture viability than pumping methods. Bidirectional Transfer Caps use a standard pump to transfer culture and/or media; come in a wide variety of tubing sizes. Transfer Caps come with multiple types of end fittings; quick connect, luer lock, and tube fusing. Please see our website htslabs.com for more details.

High cell death and a large amount of foam and/or cell clumping issues?

Cell death and foaming in the Optimum Growth™ Flasks is usually due to cell shearing. Adding up to 1% surfactant will reduce foaming and increase cell viability without stressing the cells.

How can you best use media from ThermoFisher such as F17 and its derivatives?

FreeStyle™ F17 Expression Medium contains lower amounts of pluronic than other comparable medium. Cells grown in this media may experience more shear stress due to the lower amount of surfactant. To avoid this, add in additional pluronic (ThermoFisher p/n 24040032). The recommended range of pluronic is 0.05 gm/L to 0.2 gm/L. Up to 1% simethicone from MilliporeSigma (p/n 59920C) can also be used. Either of these methods usually will work to reduce foaming and restore high culture viability.

What can I do if the doubling time for my cell culture is longer than expected when using the Optimum Growth™ Flasks?

This varies between cell types and strains, as well as with environmental conditions. If the doubling time for your culture is taking longer than expected or desired in the Thomson Optimum Growth™ Flasks, we recommend increasing the shake speed beyond our recommended speeds by 10 to 20 rpm. The reason for the increased doubling time is that the oxygen transfer rate maybe lower with higher fill volumes, and the increase in speed will compensate for this.

Disposable shake flasks are hard to remove from the sticky pad.
What do we do?

Some suggestions from people who find it too sticky:

1. Spray ethanol on the sticky pad until you reach the desired stickiness. Ethanol will lower the bonding strength, as will any alcohol.
2. Some people use rug gripper pads on top of the sticky tape.

Which transfection reagent works best with CHO & HEK293 cells?

We see that there are three classes of transfection reagents that have varying efficiencies:

Class	Example	Efficiency of Transfection
Polymers	PEI	< 65%
Cationic Liposomes	Lipofectamine™	70-95%
Electroporation	Maxcyte®	> 96%

Polymers/PEI: The most common transfection reagent used in the market. It is inexpensive, but may not lead to as high of a transfection rate, and requires higher DNA quantities. Commonly used for all small and large scale transfections.

Cationic Liposomes/Lipofectamine™: This class of transfection reagents is highly efficient and is commonly used in CHO-S, 293F and other high titer systems. Cationic Liposomes work well with our flasks. We have seen consistent transfection with great viability from TransIT-Pro®. This used with CHOgro(r) from Mirus in 24 well plates (2.5mL), small scale flasks, and production in Optimum Growth Flasks (50mL-2.5L) have given scalable usable results.

Electroporation: Most often used method for large scale, >1L transfection. Unfortunately, electroporation is not as useful for multiple transfections at one time. Customer feedback shows that stabilizing the cells with a 1%-1.2% addition of surfactant (pluronic/PF68) 30 minutes after transfection leads to higher titers and viability.



One shaker many flask sizes: consistent shake speeds across all sizes of Optimum Growth™ Flasks




Optimum Growth™ 250mL Flasks with CHO & HEK293 Cells in shaker

Thomson Instrument Company is not affiliated with Corning Life Sciences®, ThermoFisher Pluronic, MilliporeSigma, Eppendorf®, INFORS HT®, Kuhner®, Fisher, Scientific®, VWR®, Mirus or their products

Corning® vs Optimum Growth™ Flasks 2 Membrane Proteins Expressed

Data provided by Genentech part of the Roche Group



Corning® – 500mL flask, 200mL culture

Thomson – 250mL flask, 150mL culture


4mL samples purified over Ni NTA

Protein A – Membrane protein of moderate expression, 34kDa

Protein B – Membrane protein of low expression, 45kDa



12µL of elution resolved on a coomassie gel

Conclusion: Thomson flasks work at least as good as Corning® standard. Improved working volume / flask volume ratio when using Thomson.



Same Footprint-Double Volume

Optimum Growth™ 5L (3L Media) vs
Nalgene® Nunc 2.8L (1.5L Media)



Thomson Instrument Company is not affiliated with Corning Life Sciences®, Nalgene Nunc®, Genentech, NYSGRG or their products.

Regular Shake Flask vs Thomson Shake Flask

Lake Pharma

Background Information

The purpose of this study was to compare the Regular shape shake flask with Thomson shake flask for cell growth and impact on protein yield.


Conditions

HEK-293 Transient production– seeded at 2.5×10^6 cells/mL and incubated at 37°C, 5% CO₂. Supplemented with feed on day 1. Cells were harvested on day 7.




CHO Stable pool production– seeded at 0.8×10^6 cells/mL and incubated at 37°C, 5% CO₂. Supplemented with feed and glucose as needed. Temperature shift to 32°C on day 5. Cells were harvested on day 14.

Double the Volume with Thomson Optimum Growth™ Flasks

Protein Titer






293 transient production results: Case 1




293 transient production results: Case 2

900mL(WV)/1.6L Thomson vs 200mL(WV)/Regular SF1L



CHO stable pool production results



Conclusion

Transient protein expression (293)

- Both type of flasks support comparable peak VCD
- Thomson Shake Flask maintained higher VCD and viability
- Thomson Shake Flask boost protein yield more than 20%

Stable pool protein expression (CHO)

- Thomson Shake Flask supported higher peak VCD (11.4 vs 8.7×10^6 cell /mL)
- Thomson Shake Flask maintained higher culture VCD and viability
- Extended culture life for 2 days if based on the harvesting criteria
- Thomson Shake Flask results 50% titer increase



Optimum Growth™ 125mL Flasks in a shaker for cell line development



Optimum Growth™ 250mL Flasks with CHO & HEK293 Cells in shaker

Thomson Instrument Company is not affiliated with Lake Pharma, Infors or their products.


High & Low Expressing Proteins in HEK293 Cells




Low Expressing Gel

This gel shows equal bands from 5 replicates of a low expressing protein, producing roughly 10 to 20 mg/L.

1. Benchmark Pre-Stained Protein Ladder
2. Purified protein, 100 ng control
3. Purified protein, 200ng control
4. Untransfected cells, -ve control
5. +ve control
6. +ve control
7. Protein of interest, 5L Combined Flasks #1-5
8. Protein of interest, 5L Flask #1
9. Protein of interest, 5L Flask #2
10. Protein of interest, 5L Flask #3
11. Protein of interest, 5L Flask #4
12. Protein of interest, 5L Flask #5




4-20% SDS-PAGE Quick Blue Stain Commassie Gel
Expected MW of dimer 24.5 kDa Estimated expression level ~10-20 mg/L



High Expressing Gel

Thomson 5L flasks are consistently able to maximize production of your best expressers. This gel shows equal bands from 3 replicates of a high expressing protein, producing approximately 300 mg/L.

1. Benchmark Pre-Stained Protein Ladder
2. Purified mAb 100 ng control
3. Purified mAb 250 ng control
4. Purified mAb 500 ng control
5. Purified mAb 1000 ng control
6. +ve control
7. Protein of interest, 5L Flask #1
8. Protein of interest, 5L Flask #2
9. Protein of interest, 5L Flask #3
10. Protein of interest, 5L Combined Flasks #1-3



4-20% SDS-PAGE Quick Blue Stain Commassie Gel
Expected MW of dimer 159.4 kDa Estimated expression level ~300 mg/L

Conclusion

Thomson Optimum Growth™ Flasks not only ensure consistent expression from Hek293 strains, they can also increase shaker capacity. With the same footprint as a typical Corning® 3L flask and a culture volume of up to 3L, the Optimum Growth™ 5L Flask may increase production 200%, if not more, in the same space (this is construct dependent).

Most constructs express at higher levels in the Optimum Growth™ 5L flasks. This makes one Optimum Growth™ 5L equivalent to, if not greater than, two 3L flasks.

Thomson Instrument Company is not affiliated with Corning Life Sciences or their products



Key Features

- Aseptic sampling in the hood or on a shaker
- 1-way valve for simple sterile sampling
- Prevent potential contamination – no need to open the flasks

Sampling Flask

The Optimum Growth™ Flasks provide several options for aseptic sampling in all flask sizes. The 125mL, 250mL, and 500mL Optimum Growth™ Flasks have an optional 1-way valve in the vented sample cap. The 5L Optimum Growth™ Flasks has an optional 1-way valve in the side of the flask.

The aseptic sampling option is available for the 125mL, 250mL, 500mL & 5L Optimum Growth™ Flasks.

Vented Sampling Cap

The Optimum Growth™ Sampling Cap is ventilated for optimal cultivation and has a 1-way valve built into it. This reduces contamination and allows for aseptic sampling of cells while the flasks remain in the shaker, eliminating the need for transfer to the hood from the shaker.

The Optimum Growth™ Sampling Flasks come with a 0.2µm PTFE vented cap for optimum aeration during cultivation. The Thomson vented caps create a safe aseptic barrier from harmful contaminants, while the large surface area creates an optimum air exchange for cell growth.



Sample Flasks have ventilated caps for optimum aeration during cultivation.

Comparison of Old vs Thomson Sampling Method



OLD Sampling Method:

1. Remove flask from shaker
2. Spray down flask before putting in the hood
3. Placing flask in the hood
4. Removing cap
5. Take sample
6. Replace cap
7. Put back in shaker


Thomson Improved Sampling Method:





1. Connect syringe to luer lock on cap
2. Sample flask while in the shaker
3. Remove syringe

Simpler Sampling in the Shaker

Thomson's improved sampling method allows you to sample directly in the shaker without the need to remove caps or use a deteriorating needle septum.

1. Open shaker
2. Attach syringe to sampling port and withdraw the sample
3. Remove syringe & close shaker



Sampling Flask w/ Vented Cap Specifications				
Flask Size	125mL	250mL	500mL	5L
Part #	931110-SP	931111-SP	931112-SP	931116-PORT-E
Image				
Description	Optimum Growth™ 125mL Flask w/ Sampling Port Sterile	Optimum Growth™ 250mL Flask w/ Sampling Port Sterile	Optimum Growth™ 500mL Flask w/ Sampling Port Sterile	Optimum Growth™ 5L Flask w/ Sampling Port Sterile
Flask Connection	Male Luer Lock	Male Luer Lock	Male Luer Lock	Male Luer Lock
Tubing Connection	n/a	n/a	n/a	n/a
Sample Tubing Volume	163µL	218µL	313µL	381µL
Top Style	Threaded	Threaded	Threaded	Threaded
Cap Material	PP (polypropylene)	PP (polypropylene)	PP (polypropylene)	PP (polypropylene)
Sterility	10 ⁻⁶	10 ⁻⁶	10 ⁻⁶	10 ⁻⁶
Air Filter/Ventilation	0.2µm PTFE for cultivation & pressure relief	0.2µm PTFE for cultivation & pressure relief	0.2µm PTFE for cultivation & pressure relief	0.2µm PTFE for cultivation & pressure relief
Qty/Case	50	50	25	4



Key Features

- Sterile tube fuse inoculation
- 1-way valve sampling allowing for a simple use in a shaker
- Eliminate the need for under 20L cell bags
- Replaces current process requiring tube fusing for inoculation
- Consistently allows for multiple day additions
- Allows for simple use within a shaker
- Great for aseptic manufacturing; never any opening needed preventing potential contamination

Multiported Optimum Growth™ Transfer & Feed Flask

The multiported flask is completely aseptic, making it the perfect start for initial seed cultures that inoculate bioreactors in multiple stages of clinical drug production. Other uses for the Multiported Optimum Growth™ Flasks include maintaining processes requiring long term duration cell culture and other manufacturing functions.

Born out of necessity, the Multiported Optimum Growth™ Flask has replaced the need for bags to start seed cultures for the inoculation of bioreactors in biopharmaceutical companies




5L TRT Flask, Transfer/Feed side





5L TRT Flask, Sample port side

Multiported Optimum Growth™ Transfer & Feed Explained

One way sample port allows safe sampling of cells while the flasks remain in the shaker.



Two way port allows feeding of media during growth and transfer to bioreactor for seeding a larger production growth.

Multiported Optimum Growth™ Transfer & Feed Flask Specifications		
Flask Size	1.6L	5L
Part #	931113-PORT-TRT	931116-PORT-TRT-F
Image		
Description	Multiported Optimum Growth™ 1.6L Flask w/ Bidirectional Transfer/Feed & Sampling Ability Sterile	Multiported Optimum Growth™ 5L Flask w/ Bidirectional Transfer/Feed & Sampling Ability Sterile
Flask Connection	Male Luer Lock	Male Luer Lock
Tubing Connection	Tube Fuse	Tube Fuse
Tubing Diameter	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)
Tubing	Chemically resistant, heat sealable, flexible	Chemically resistant, heat sealable, flexible
Tubing Length	24" (609.6mm)	24" (609.6mm)
Sample Tubing Volume	326µL	381µL
Top Style	Threaded	Threaded
Cap Material	PP (polypropylene)	PP (polypropylene)
Sterility	10 ⁻⁶	10 ⁻⁶
Air Filter Ventilation	0.2µm PTFE vent filter	0.2µm PTFE vent filter
Qty/Case	12	4



Key Features

- Transfer from an Optimum Growth™ Flask to another Flask, cell bag or bioreactor
- Gravity feed keeps cells stress free
- Allows for simple aseptic transfer outside of hood

Inversion Transfer Cap For 1.6L – 5L

The Sterile Inversion Transfer Caps (patented), allow for the efficient aseptic transfer of media or cells to cell bags, bioreactors, or flasks (from all manufacturers).

The Inversion Transfer Cap works with the 1.6L, 2.8L and 5L Thomson Optimum Growth™ Flasks (patented). Simply, replace the culture cap with the Transfer Cap and connect to your vessel of choice.

How to Transfer?

To transfer, invert the Optimum Growth™ Flask and let gravity do the rest. The Optimum Growth™ Flask and Transfer Cap System eliminates the need to transfer cells to an intermediate flask for scale up. Inversion or gravity feed has the lowest shear force of any transfer method available.

The Inversion Transfer Cap in conjunction with the 1.6L, 2.8L and 5L Thomson Optimum Growth™ Flasks (patented) product line can be used for reagent addition/seeding of larger bioreactors cell bags, or other liquid media transfers.


Connection Options




Quick Connect




Luer Lock



Quick Connect



Tube Fuse




Transfection/Scale-Up from Optimum Growth™ 5L Flask to Bag



Transfer from Flask to Bag









Completed Aseptic Transfer outside of hood



Stand with Ring are available for hands free transfer.

Inversion Transfer Cap Accessories				
Flask Size	1.6L	1.6L	5L	5L
Part #	931609	931700	931606	931607
Image				
Description	Stand with Ring for Inverting Optimum Growth™ 1.6 & 2.8L Flasks with Inversion Transfer Cap	Ring Only for Inverting Optimum Growth™ 1.6 & 2.8L Flasks requires a 22" stand	Stand with Ring for Inverting Optimum Growth™ 5L Flasks with Inversion Transfer Cap	Ring Only for Inverting Optimum Growth™ 5L Flask requires a 22" stand
Stand Height	22"	n/a	22"	n/a
Ring Diameter	5"	5"	7"	7"
Qty/Case	1	1	1	1

Inversion Transfer Cap Specifications				
Flask Size	1.6L & 2.8L	1.6L & 2.8L	1.6L & 2.8L	1.6L & 2.8L
Part #	931706	931710	931705	931708
Image				
Description	Inversion Transfer Cap for Optimum Growth™ 1.6 & 2.8L Flask	Inversion Transfer Cap for Optimum Growth™ 1.6 & 2.8L Flask	Inversion Transfer Cap for Optimum Growth™ 1.6 & 2.8L Flask	Inversion Transfer Cap for Optimum Growth™ 1.6 & 2.8L Flask
Tubing Connection	7/16" (11.1mm) Male Connection	Tube Fuse/Female Luer Lock	Tube Fuse	Tube Fuse
Tubing Diameter	n/a	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)	C-Flex® 24 ID: 3/16" (4.76mm), OD: 7/16" (11.1mm)
Tubing	n/a	Chemically resistant, heat sealable, flexible	Chemically resistant, heat sealable, flexible	Chemically resistant, heat sealable, flexible
Tubing Length	n/a	24" (609.6mm)	24" (609.6mm)	24" (609.6mm)
Top Style	Threaded	Threaded	Threaded	Threaded
Cap Material	PP (polypropylene)	PP (polypropylene)	PP (polypropylene)	PP (polypropylene)
Sterility	10 ⁻⁶	10 ⁻⁶	10 ⁻⁶	10 ⁻⁶
Air Filter Ventilation	0.2µm PTFE vent filter	0.2µm PTFE vent filter	0.2µm PTFE vent filter	0.2µm PTFE vent filter
Qty/Case	4	4	4	4

Inversion Transfer Cap Specifications					
Flask Size	5L	5L	5L	5L	5L
Part #	931594	931596	931616	931595	931598
Image					
Description	Inversion Transfer Cap for Optimum Growth™ 5L Flask	Inversion Transfer Cap for Optimum Growth™ 5L Flask	Inversion Transfer Cap for Optimum Growth™ 5L Flask	Inversion Transfer Cap for Optimum Growth™ 5L Flask	Inversion Transfer Cap for Optimum Growth™ 5L Flask
Tubing Connection	1/4" (6.35mm) Barb	7/16" (11.1mm) Barb	Tube Fuse/Female Luer Lock	Tube Fuse	Tube Fuse
Tubing Diameter	n/a	n/a	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)	C-Flex® 24 ID: 3/16" (4.76mm), OD: 7/16" (11.1mm)
Tubing	n/a	n/a	Chemically resistant, heat sealable, flexible	Chemically resistant, heat sealable, flexible	Chemically resistant, heat sealable, flexible
Tubing Length	n/a	24" (609.6mm)	24" (609.6mm)	24" (609.6mm)	24" (609.6mm)
Top Style	Threaded	Threaded	Threaded	Threaded	Threaded
Cap Material	PP (polypropylene)	PP (polypropylene)	PP (polypropylene)	PP (polypropylene)	PP (polypropylene)
Sterility	10 ⁻⁶	10 ⁻⁶	10 ⁻⁶	10 ⁻⁶	10 ⁻⁶
Air Filter Ventilation	0.2µm PTFE vent filter	0.2µm PTFE vent filter	0.2µm PTFE vent filter	0.2µm PTFE vent filter	0.2µm PTFE vent filter
Qty/Case	4	4	4	4	4



Key Features

- Tube fuse for aseptic transfer of media or cells into or out of flasks, cell bags, & bioreactors
- Bidirectional Transfer Cap with a downstem allows for pumping of liquid into or out of an Optimum Growth™ Flask
- Available for 1.6L & 5L Optimum Growth™ Flasks

Bidirectional Transfer Cap For 1.6-L – 5L

The Sterile Optimum Growth™ Bidirectional Transfer Caps (patented) with downstem, allow for easy aseptic transfer of media or cells into and/or out of cell bags, bioreactors, and flasks (from all manufacturers).

The Bidirectional Transfer Cap with downstem works with the 1.6L, 2.8L and 5L Thomson Optimum Growth™ Flasks (patented) and a peristaltic pump. Simply, replace the culture cap with the Transfer Cap and connect to your vessel of choice.


How to Transfer or Feed

1. Replace with Bidirectional Transfer Cap
2. Connect to your receiving vessel by Tube Fusing or using our Luer Lock option
3. Place the silicone tubing in the peristaltic pump head
4. Liquid can then be pumped either into or out of the flask


Connection Options



Luer Lock




Tube Fuse



Optimum Growth™ Flasks (patented) in conjunction with the Transfer Cap System eliminates the need to move cells to an intermediate transfer for scale up or seed cultures. The ability to pump into the Optimum Growth™ Flask (patented) makes filling with media from a bulk source a simple aseptic method. The Transfer Cap in conjunction with the Thomson Optimum Growth™ Flask product line can be used for reagent addition, seeding of larger bioreactors or cell bags, pumping of media into flasks from large drums or bags of media, and other liquid media transfers into and out of bioreactors.

The Optimum Growth™ Flasks come in multiple sizes of 125mL, 250mL, 500mL, 1.6L, 2.8L and our popular 5L.



Transfer from flask to bioreactor with pump

Thomson Instrument Company is not affiliated with New Brunswick

Bidirectional Transfer Cap Specifications					
Flask Size	1.6L	1.6L	2.8L	5L	5L
Part #	931702	931704	931804	931618	931614
Image					
Description	Bidirectional Transfer Cap for Optimum Growth™ 1.6L Flask	Bidirectional Transfer Cap for Optimum Growth™ 1.6L Flask	Bidirectional Transfer Cap for Optimum Growth™ 2.8L Flask	Bidirectional Transfer Cap for Optimum Growth™ 5L Flask	Bidirectional Transfer Cap for Optimum Growth™ 5L Flask
Tubing Connection	Tube Fuse/Luer Lock	Tube Fuse	Tube Fuse	Tube Fuse/Luer Lock	Tube Fuse
Tubing Diameter	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)	C-Flex® 16 ID: 1/8" (3.1mm), OD: 1/4" (6.35mm)
Tubing	Chemically resistant, heat sealable, flexible	Chemically resistant, heat sealable, flexible	Chemically resistant, heat sealable, flexible	Chemically resistant, heat sealable, flexible	Chemically resistant, heat sealable, flexible
Tubing Length	24" (609.6mm)	24" (609.6mm)	24" (609.6mm)	24" (609.6mm)	24" (609.6mm)
Top Style	Threaded	Threaded	Threaded	Threaded	Threaded
Cap Material	PP (polypropylene)	PP (polypropylene)	PP (polypropylene)	PP (polypropylene)	PP (polypropylene)
Sterility	10 ⁻⁶	10 ⁻⁶	10 ⁻⁶	10 ⁻⁶	10 ⁻⁶
Air Filter Ventilation	0.2µm PTFE vent filter	0.2µm PTFE vent filter	0.2µm PTFE vent filter	0.2µm PTFE vent filter	0.2µm PTFE vent filter
Qty/Case	8	8	8	8	8



Key Features

- 0.2µm filter, eliminates multiple filtration steps and in some cases centrifugation*
- Significant time savings versus traditional spin down technique
- Cell clarification of low or high density cultures of CHO stable, CHO transient, HEK293, hybridoma, and other mammalian cell lines
- Eliminates transfer steps: The Rapid Clear® Cap screws directly onto the Optimum Growth™ Flask
- Secondary cap attaches to a new Optimum Growth™ Flask or to a storage container with a Luer lock




Current Method




New Method

Rapid Clear®

The Thomson Optimum Growth™ Family of products is expanding into downstream processing with a revolutionary new technology that allows high speed clarification of cellular material. The filters work by simply switching the vented Optimum Growth™ Cap to the Rapid Clear® Cap. The only equipment required is a peristaltic pump. In initial testing the units can filter up to 2.5L-3L of high density culture, 6×10^6 - 20×10^6 @ >70% viability, in approximately 18 minutes. This technique will transform the time consuming and laborious task of harvesting cells to a rapid walk away procedure.





Rapid Clear® Cap 3000 Specifications


Flask Size	2.8L & 5L
Part #	788116
Image	
Description	High-speed clarification of 2L–4L of cell culture Sterile
Tubing Connection	Transfer Cap
Tubing Diameter	Size 15 silicone tubing, ID: 3/16" (4.76mm), OD: 7/16" (11.1mm)
Tubing	Chemically resistant, flexible
Tubing Length	48" (1219.2mm)
Cap Material	PP (polypropylene)
Qty/Case	8

Part Numbers & Specifications





Optimum Growth™ Flask Specifications						
Flask Size	125mL	250mL	500mL	1.6L	2.8L	5L
Part #	931110	931111	931112	931113	931114	931116
Image						
Description	Optimum Growth™ 125mL Flask w/ Vent Cap Sterile	Optimum Growth™ 250mL Flask w/ Vent Cap Sterile	Optimum Growth™ 500mL Flask w/ Vent Cap Sterile	Optimum Growth™ 1.6L Flask w/ Vent Cap Sterile	Optimum Growth™ 2.8L Flask w/ Vent Cap Sterile	Optimum Growth™ 5L Flask w/ Vent Cap Sterile
Qty/Case	50	50	25	12	6	4





Sampling Flask w/ Vented Cap Specifications				
Flask Size	125mL	250mL	500mL	5L
Part #	931110-SP	931111-SP	931112-SP	931116-PORT-E
Image				
Description	Optimum Growth™ 125mL Flask w/ Sampling Port & Vent Cap Sterile	Optimum Growth™ 250mL Flask, w/ Sampling Port & Vent Cap Sterile	Optimum Growth™ 500mL Flask, w/ Sampling Port & Vent Cap Sterile	Optimum Growth™ 5L Flask, w/ Sampling Port & Vent Cap Sterile
Qty/Case	50	50	25	4





Multiported Flask w/ Vented Cap Specifications	
Flask Size	1.6L
Part #	931113-PORT-TRT
Image	
Description	Multiported Optimum Growth™ 1.6L Flask w/ Bidirectional Transfer/Feed & Sampling Ability Sterile
Qty/Case	12
Flask Size	5L
Part #	931116-PORT-E-TRT
Image	
Description	Multiported Optimum Growth™ 5L Flask w/ Bidirectional Transfer/Feed & Sampling Ability Sterile
Qty/Case	4

Rapid Clear® Cap 3000 Specifications	
Flask Size	2.8L & 5L
Part #	788116
Image	
Description	High-speed clarification of 2–4L of cell culture for use with 2.8 & 5L Optimum Growth™ Flasks
Qty/Case	8



Part Numbers & Specifications

Inversion Transfer Cap Specifications				
Flask Size	1.6L & 2.8L	1.6L & 2.8L	1.6L & 2.8L	1.6L & 2.8L
Part #	931706	931710	931705	931708
Image				
Description	Inversion Transfer Cap for Optimum Growth™ 1.6 & 2.8L Flask, 7/16" Male Connection Sterile	Inversion Transfer Cap for Optimum Growth™ 1.6 & 2.8L Flask, with 2' C-Flex® 16 Tubing w/ Luer Lock Sterile	Inversion Transfer Cap for Optimum Growth™ 1.6 & 2.8L Flask, 2' C-Flex® 16 (1/4") Tubing Sterile	Inversion Transfer Cap for Optimum Growth™ 1.6 & 2.8L Flask, 2' C-Flex® 24 Tubing Sterile
Qty/Case	4	4	4	4





Inversion Transfer Cap Specifications				
Flask Size	5L	5L	5L	5L
Part #	931594	931596	931616	931595
Image				
Description	Inversion Transfer Cap for Optimum Growth™ 5L Flask, 1/4" Barb Connection Sterile	Inversion Transfer Cap for Optimum Growth™ 5L Flask, 7/16" Male Connection Sterile	Inversion Transfer Cap for Optimum Growth™ 5L Flask, w/ 2' C-Flex® 16 (1/4") Tubing w/ Luer Lock Sterile	Inversion Transfer Cap for Optimum Growth™ 5L Flask, 2' C-Flex® 16 (1/4") Tubing Sterile
Qty/Case	4	4	4	4

Bidirectional Transfer Cap Specifications				
Flask Size	1.6L	1.6L	2.8L	5L
Part #	931702	931704	931804	931618
Image				
Description	Bidirectional Transfer Cap for Optimum Growth™ 1.6L Flask, 2' Tubing w/ Luer Lock Sterile	Bidirectional Transfer Cap for Optimum Growth™ 1.6L Flask, 2' C-Flex® 16 (1/4") tubing Sterile	Bidirectional Transfer Cap for Optimum Growth™ 2.8L Flask, 2' C-Flex® 16 (1/4") tubing Sterile	Bidirectional Transfer Cap for Optimum Growth™ 5L Flask, 2' Tubing w/ Luer Lock Sterile
Qty/Case	8	8	8	8

Accessories

Optimum Growth™ Accessories		
Flask Size	125mL	250mL
Part #	1212900	1212905
Image		
Description	8 Position Carrier for Optimum Growth™ 125mL Flask (931110)	8 Position Carrier for Optimum Growth™ 250mL Flask (931111)
Flask Capacity	8	8
Dimensions	10.75" x 5"	13.4" x 6"
Qty/Case	1	1

Optimum Growth™ Accessories					
Flask Size	125mL	250mL	500mL	1.6L, 2.8L & 5L	1.6L, 2.8L & 5L
Part #	899110	899111	899112	899116	899600-B
Image					
Description	Thomson Vent Caps for Optimum Growth™ 125mL Flasks Sterile	Thomson Vent Caps for Optimum Growth™ 250mL Flasks Sterile	Thomson Vent Caps for Optimum Growth™ 500mL Flasks Sterile	Thomson Vent Caps for Optimum Growth™ 1.6, 2.8 & 5L Flasks Sterile	Thomson Solid Cap for Optimum Growth™ 1.6, 2.8 & 5L Flasks Sterile
Qty/Case	50	50	25	24	24

Inversion Transfer Cap Accessories				
Flask Size	1.6L	1.6L	5L	5L
Part #	931609	931700	931606	931607
Image				
Description	Stand with Ring for Inverting Optimum Growth™ 1.6 & 2.8L Flasks with Inversion Transfer Cap	Ring Only for Inverting Optimum Growth™ 1.6 & 2.8L Flasks requires a 22" stand	Stand with Ring for Inverting Optimum Growth™ 5L Flasks with Inversion Transfer Cap	Ring Only for Inverting Optimum Growth™ 5L Flask requires a 22" stand
Qty/Case	1	1	1	1



THOMSON
MADE IN USA

